

VisionArray® HPV Chip 1.0

Results of the <u>European HPV DNA Test External</u> <u>Quality Assurance Scheme (EHEQAS)</u>



Mendel Center for Biomedical Sciences • www.mendelcenter.org/eheqas

Dear Readers,

With this new issue of our **ZYTONEWS** we would like to summarize the results of the **European HPV DNA Test External Quality Assurance Scheme (EHEQAS)** results achieved in April 2016 and May 2017. EHEQAS is an EQA scheme initiated in 2003 and formally founded in 2006. It already includes experts and diagnostic laboratories from Croatia, Cyprus, Czech Republic, Germany, Greece, Hungary, Netherlands, and Portugal.

Enjoy reading, Yours ZYTONEWS TEAM



Aim & Method of the Quality Testing

EHEQAS is an **External Quality Assurance Scheme** to assess and improve the quality of HPV DNA testing laboratories.

Any European laboratory performing HPV tests may participate. After registration the participants receive 5-7 samples for HPV detection 1-2 times per year. It is considered that participants should be able to report results to the coordinator within 3 weeks from receiving the samples. Subsequently, test reports of the results are evaluated by the coordinating team and the results and methods obtained from all participants are compared. Consensus results are issued and announced to the participants. Successful participation in EHEQAS is extremely helpful for the laboratories to verify and validate the results produced.

Overview of applied HPV Detection Methods

Test Assay	Method	Genotyping
Chipron HPV 3.5	PCR/reverse hybridization on LCD-Chip	Ø
CLART® HPV	PCR/hybridization on a microarray surface, followed by an enzymatic colour reaction	Ø
Digene HC2 HPV DNA Test	DNA-RNA hybrids are captured by antibodies, AP-detection with chemiluminescent substrate	-
INNO-LiPA® HPV Genotyping Extra II	PCR/reverse hybridization on membrane strips, followed by an enzymatic colour reaction	Ø
VisionArray [®] HPV Chip 1.0	PCR/hybridization on the VisionArray® HPV Ch followed by an enzymatic colour reaction	ip, Ø

Assessment of the Results

After receiving the test reports, the coordinating team is evaluating the results for each sample. Every tested sample can be assessed with a maximum of 10 marks:

- 5 of these 10 marks are awarded for a detection of the correct presence or absence of HPV DNA (Possible results: Negative, positive, insufficient sample, PCR inhibitor).
- The other 5 marks are awarded for a correct HPV genotyping. The correct HPV genotype as well as the detection of a possible multi-infection is of importance.

Quality Control Trial I - April 2016

The results of 8 laboratories participating in EHEQAS April 2016 are summarized. Each laboratory received the same 5 or 7 EHEQAS samples (same requirements for all laboratories) which should be processed by applying their routine HPV testing method. If all 5 and 7 samples were evaluated successfully, the maximum score achievable was 50 and 70, respectively. Two laboratories achieved the maximum score: Genomica HPV 3 and **VisionArray**[®] **HPV Chip 1.0** as HPV testing method.



Results: April 2016

Conclusion

ZytoVision VisionArray[®] **HPV Chip 1.0** was used as testing method in 1 laboratory. ZytoVision achieved a full score of 50 marks (100%).

Quality Control Trial II - May 2017

The results of 10 laboratories participating in EHEQAS May 2017 are summarized. Each laboratory received the same 5 or 7 EHEQAS samples (same requirements for all laboratories) which should be processed by applying their routine HPV testing method. If all 5 and 7 samples were evaluated successfully, the maximum score achievable was 50 and 70, respectively. Two laboratories achieved the maximum score: InnoLipa HPV genotyping extra strip assay and

VisionArray[®] HPV Chip 1.0 as HPV testing method.



Results: May 2017

Conclusion

The **VisionArray**[®] **HPV Chip 1.0** was used as testing method in 5 laboratories. ZytoVision achieved a full score of 50 marks (100%). In 3 other laboratories using **VisionArray**[®] for testing, the obtained results were assessed as excellent (96-100%) demonstrating the high capability and reproducibility of the system.

Summary of Conclusions

- The **ZytoVision labaratory** participated successfully in EHEQAS in April 2016 and May 2017 with a **100% assessment**.
- Other laboratories started to use VisionArray[®] as testing method with excellent EHEQAS results (96-100%).
- The VisionArray[®] system reached 100% grading in both EHEQAS HPV quality control trials, demonstrating the high reliability, sensitivity and reproducibility of this HPV test.

VisionArray® – Method Description



Product Information

Simultaneous genotyping of 41 different HPV types - all certified for in vitro diagnostic use

- All capture sequences and positive controls are set up on the VisionArray® HPV Chips as duplicates
- High sensitivity and specificity
- Quick & easy 1 hour protocol
- Automated evaluation using the VisionArray[®] Analyzer Software simple visualization & quick analysis in just a few minutes



(Prod. No.	Product	Tests	_
	VA-0001-10	VisionArray HPV Chip 1.0 Ind. 10 pieces CE IVD	10	
	VA-0001-50	VisionArray HPV Chip 1.0 Ind. 5x 10 pieces CE IVD	50	
	VA-0002-10	Vision <i>Array</i> HPV High Risk Chip 1.0 Ind. 10 pieces C€ IVD	10	
	VA-0002-50	VisionArray HPV High Risk Chip 1.0 Incl. 5x 10 pieces CE IVD	50	

Product Information

VisionArray[®] Analyzer Software

- Simple visualization and quick analysis of the Vision Array® HPV Chip data
- Analysis of the VisionArray® HPV Chip and the report of the results can be achieved in just a few minutes
- Program navigation is very easy and intuitive
- Scans are stored in an **integrated database** on the enclosed external hard drive including all sample and chip data



Prod. No. Product

E-4060-1 Vision Array Analysis Package C E IVD Ind. Scanner 8100; Side Holder; Hand Scanner; PC with preinstalled Vision Array Analyzer Software; USB-Hub; External Hard Drive; Computer Mouse

VisionArray® Detection Kit

For hybridization and detection of PCR products on VisionArray® Chips

Prod. No.	Product	Tests
VK-0003-50	VisionArray Detection Kit CE IVD	50
	Incl. Hybridization Solution, 1 ml; Detection Solution, 5 ml; Blue Spot Solution, 5 ml; 100x Wash Buffer, 250 ml	
E-4051-1	Mini Slide Centrifuge	
	Incl. 2 place slide rotor, two slide holders	

VisionArray[®] DNA Extraction Kits

For isolation of genomic DNA from FFPE as well as liquid based cytology specimens

$\left(\right)$	Prod. No.	Product	Tests	
	VI-0001-50	Vision <i>Array</i> FFPE DNA Extraction Kit Incl. Paraffin Dissolver; Tissue Lysis Buffer, Decrosslink Buffer; DNA Wash Buffer, Proteinase K. Proteinase K Buffer, Elution Buffer, Columns; Collection Tubes	50	
	VI-0002-50	Vision <i>Array</i> Cytology DNA Extraction Kit Incl. Pre-Lysis Buffer; Cell Lysis Buffer, DNA Wash Buffer, Proteinase K; Proteinase K Buffer, Elution Buffer, Columns; Collection Tubes	50	

VisionArray[®] PCR Reagents

For contamination-free amplification and biotinylation of target sequences with a high quality heat stable Taq polymerase

Prod. No.	Product	Tests	
VP-0001-50	Vision <i>Array</i> HPV Primer Kit CE IVD Ind. HPV Primer Mix; dNTP/dUTP Solution	50	
VE-0001-100	Vision <i>Array</i> PreCise Taq DNA Polymerase CE IVD Incl. VisionArray PreCise Taq DNA Polymerase; Reaction Buffer, 10x; PreCise MgCL ₂ , 25 mM	100	
VE-0002-100	Vision <i>Array</i> Uracil-DNA Glycosylase CE IVD Incl. VisionArray Uracil-DNA Glycosylase (UDG); UDG Reaction Buffer, 10x	100	
ES-0007-50	Vision <i>Array</i> HPV PreCise Master Mix CE IVD	50	



Molecular diagnostics simplified

ZytoVision GmbH · Fischkai 1 27572 Bremerhaven · Germany Phone: +49 (0)471/4832-300 Fax: +49 (0)471/4832-509 info@zytovision.com

For further information, please contact us.

