



VisionArray®

Arrays for DNA analysis

Genotyping of High Risk Human Papillomavirus

VisionArray® High Risk HPV

Your Complete
Solution from
DNA Extraction
to Analysis!

Simultaneous Detection of 24 HPV Types

High Sensitivity – Easy to Handle

ZYTOVISION
Molecular diagnostics simplified

VisionArray® Arrays for DNA analysis

VisionArray® HPV High Risk Chip 1.0



Prod. No.	Product	Tests
VA-0002-10	VisionArray HPV High Risk Chip 1.0 CE IVD Incl. 10 pieces	10
VA-0002-50	VisionArray HPV High Risk Chip 1.0 CE IVD Incl. 5x 10 pieces	50

VisionArray® DNA Extraction, PCR, and Detection



Prod. No.	Product	Tests
VI-0001-50	VisionArray FFPE DNA Extraction Kit Incl. Paraffin Dissolver; Tissue Lysis Buffer; Decrosslink Buffer; DNA Wash Buffer; Proteinase K; Proteinase K Buffer; Elution Buffer; Columns; Collection Tubes	50
VI-0002-50	VisionArray Cytology DNA Extraction Kit Incl. Pre-Lysis Buffer; Cell Lysis Buffer; DNA Wash Buffer; Proteinase K; Proteinase K Buffer; Elution Buffer; Columns; Collection Tubes	50



Prod. No.	Product	Tests
VP-0001-50	VisionArray HPV Primer Kit 2.0 CE IVD Incl. HPV Primer Mix 2.0; dNTP/dUTP Solution	50
VE-0001-100	VisionArray PreCise Taq DNA Polymerase CE IVD Incl. VisionArray PreCise Taq DNA Polymerase; PreCise Reaction Buffer, 10x; PreCise MgCl ₂ , 25 mM	100
VE-0002-100	VisionArray Uracil-DNA Glycosylase CE IVD	100
ES-0007-50	VisionArray HPV PreCise Master Mix CE IVD Containing HPV Primer Mix 2.0; dNTP/dUTP Solution; VisionArray PreCise Taq DNA Polymerase; PreCise Reaction Buffer, 10x; PreCise MgCl ₂ , 25 mM; VisionArray Uracil-DNA Glycosylase	50



Prod. No.	Product	Tests
VK-0003-50	VisionArray Detection Kit CE IVD Incl. Hybridization Solution, 1 ml; Detection Solution, 5 ml; Blue Spot Solution, 5 ml; 100x Wash Buffer, 250 ml	50
E-4051-1	Mini Slide Centrifuge	

VisionArray® Analysis Package



Prod. No.	Product
E-4060-1	VisionArray Analysis Package CE IVD Incl. Scanner 8100; Slide Holder; Hand Scanner; PC with preinstalled VisionArray Analyzer Software; USB-Hub; External Hard Drive; Computer Mouse

CE IVD only available in certain countries. All other countries research use only! Please contact your local dealer for more information.

Manufactured by: 42 life sciences GmbH & Co. KG · Fischkai 1 · 27572 Bremerhaven · Germany
Distributed by: ZytoVision GmbH · Fischkai 1 · 27572 Bremerhaven · Germany · www.zytovision.com · info@zytovision.com
VisionArray® is a trademark of 42 life sciences GmbH & Co. KG

ZYTOVISION
Molecular diagnostics simplified